

ImageJ: a useful image  
reprocessing tool

# Contents


- Remeasure distance on an SEM (JPG)
- Calculate open area
- Construct 3D images

# Remeasure an SEM image

- Download from <https://imagej.nih.gov/ij/download.html>
- Alternative versions available like Fuji.

# Remeasure an SEM image


Open an image



# Remeasure an SEM image


Use 'line' tool to  
draw over the scale

Hold shift for a  
horizontal line




# Remeasure an SEM image

Analyze -> Set Scale...



# Remeasure an SEM image

Tell ImageJ what the known scale length is



# Remeasure an SEM image

OR, even better! Just use the AutoScale marco!

It finds the scale mark, reads the text, set it up for you!

The only known issue: can sometimes tell you the unit is `Âµm` rather than `µm`...

This is because `µm` in Unicode is `Âµm` in UTF-8



So check your encoding...

```
////////// begin macro
requires( "1.51t" );
run( "8-bit" );
b = 10;
setScaleFromBar( b );
exit();
////////// end macro
////////// begin functions
function diff( aa ) {
    for ( i = 1; i < aa.length; i++ ) {
        aa[i-1] = aa[i] - aa[i-1];
    }
    aa[aa.length-1] = 0;
}
function getBeginEnd( aaa ) {
    mi = Array.findMinima( aaa, 0 );
    mx = Array.findMaxima( aaa, 0 );
    if ( mx[0] < mi[0] ) aaa[1] = mx[1]; else aaa[1] = mx[0];
    if ( mi[0] > aaa[1] ) aaa[0] = mi[1]; else aaa[0] = mi[0];
}
function setScaleFromBar( bb ) {
    run( "Set Scale...", "distance=0 known=0 pixel=1 unit=pixel" );
    ww = getWidth();
    makeRectangle( 0, bb, ww, getHeight() - 2 * bb );
    setKeyDown( "alt" );
    a = getProfile();
    setKeyDown( "none" );
    min = Array.findMinima( a, 0 );
    if ( min[0] > min[1] ) yy = min[1]; else yy = min[0];
    pel = abs( min[1] - min[0] );
    a = Array.slice( a, yy + pel );
    diff( a );
    getBeginEnd( a );
    makeRectangle( 0, yy + bb + pel + a[0], ww, abs( a[0] - a[1] ) );
    a = getProfile();
    run( "Select None" );
    diff( a );
    getBeginEnd( a );
    pel = abs( a[0] - a[1] );
    run( "Set Scale...", "distance=" + pel + " known=500 unit=µm" );
}
////////// end functions
```



# Remeasure an SEM image

Measure a gap of  
your choice





# Remeasure an SEM image

Measure a gap of  
your choice

Then Ctrl+M

Results						
File	Edit	Font	Results			
	Area	Mean	Min	Max	Angle	Length
1	147140.338	56.277	54	83	0	9253.96
2	148400.101	55.336	55	71	0	9333.33




Calculate open area

# Calculate open area

Customer told us the OA is 40%

We are not convinced...




# Calculate open area

Name the file `1.jpg`

Put it in the same folder as `OA.exe` (or `OA.py` if you have all the libs installed)


OA.exe is available my SEM folder



# Calculate open area

Run OA.exe

A window will pop up showing the processed binary image




# Calculate open area

Run OA.exe in cmd/powershell

A window will pop up showing the processed binary image

Command line will tell you `Open area is 37.69%`




# Calculate open area

Alternatively, use ImageJ

Binarise the image, denoise it by using remove outliers

Find the result in Image-Adjust-Threshold

Result is 39.03% (thx to denoise algorithm)






# Calculate open area

Marco here:


```
setOption("BlackBackground", false);  
run("Make Binary");  
run("Remove Outliers...", "radius=2 threshold=50  
which=Bright");  
run("Remove Outliers...", "radius=2 threshold=50  
which=Bright");  
setAutoThreshold("Default");  
//run("Threshold...");
```



Construct 3D images

# Construct 3D images – single SEM

This method reconstruct 3D image using the greyscale info from a single SEM image.  
Simple use plugins->3D-> Interactive surface plot



# Construct 3D images - Stack

This is yet another method to give you beautiful 3D images

First take a `stack' of images. Here im using an optical scope as an example.

I take a photo with increasing Z. This gives photos of different focusing. Ideally the Z should be equally spaced.

Use File->Import ->Image Sequence

# Construct 3D images - Stack

Use plugins->3D->Stack 3D viewer

Example taken from my thesis...

